

Chemical Resistance Technical Bulletin

OMNOVA Solutions is a leading designer and manufacturer of decorative laminates used in healthcare environments and applications, ranging from cabinetry to over-the-bed trays. Since their introduction, **surf(x)**[®] 3D Laminates have developed a strong reputation in healthcare environments. **surf(x)** 3D Laminates can be membrane-pressed or vacuum-formed to contoured surface profiles, eliminating the need for T-molding, edgebanding and visible seams. Elimination of seams helps mitigate the accumulation of potential contaminants and simplifies cleaning.

A stringent program of cleaning and disinfecting has been proven to reduce the incidence of hospital acquired infections. OMNOVA's **surf(x)** 3D Laminates can be safely cleaned and disinfected with the following cleaners and disinfectants. In laboratory testing using both the NEMA LD 3.4 and ASTM D3450 test methods, **surf(x)** 3D Laminates were shown to have no discoloration or damage to the finish when cleaned and disinfected with the following products:

- Asepticare TB-II
- Bleach 5:1
- Bleach-Rite[®] disinfecting spray
- Cavicide[®] surface disinfectant
- Citrace[®] germicide
- Clorox[®] germicidal wipes
- Dispatch[®] spray hospital cleaner disinfectant with bleach
- Fade-A-Dyne[®] stain remover
- Lysol[®] spray disinfectant
- Precise™ hospital cleaner
- PDI[®] Super Sani-Cloth[®] germicidal wipes
- PDI Sani-Cloth HB germicidal wipes
- PDI Sani-Cloth Plus germicidal wipes
- SaniZide Plus[®] germicidal solution
- Virex[®] II 256 cleaner
- Virox[®] AHP 5 disinfectant cleaner surface wipes

Trademarks are the property of their respective owners.

For the latest designs, more information or to place an order, contact OMNOVA Solutions:
866-332-5226 • laminates@omnova.com • www.omnova.com/laminates



World Headquarters
175 Ghent Road
Fairlawn, OH 44333
United States
PHONE +1 330 869-4200
FAX +1 330 869-4212